

# CREDCHROM HPLC COLUMN



**CREDENCE SCIENTIFIC**

# PREPARATIVE HPLC COLUMN

Credence Scientific supply high-performance and economical types of Preparative HPLC columns to meet your multiple needs. high-performance columns use materials all from the U.S., Japan, and Europe to guarantee your usage. All of CredChrom columns are high Performance guaranteed.



## HPLC Preparative Columns

Inner diameter: 7.8mm / 10mm / 20mm / 21.2mm / 30mm / 50 mm

Length: 50 - 500 mm

## PREP HPLC Column Introduction


Chemistry	Particle Size	Pore Size	Specific surface area	Carbon content	End capped	Usage pH
C18	5/8/10/15/20/50/75 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	16%	Yes	pH 2-8
C18H	5/8/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	20%	Yes	pH 2-11
C8M	5/8/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	12%	Yes	pH 2-8
C8X	5/8/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	9%	Yes	pH 2-11
C4	5/8/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	3%	Yes	pH 2-8
C18Bio	5/10 $\mu$ m	300 $\text{\AA}$	120m <sup>2</sup> /g	8%	Yes	pH 2-8
C8Bio	5/10 $\mu$ m	300 $\text{\AA}$	120m <sup>2</sup> /g	5%	Yes	pH 2-8
Diol	5/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	8%	Yes	pH 2-8
NH2	5/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	5%	Yes	pH 2-8
CN	5/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	7%	Yes	pH 2-8
Phenyl	5/10 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	8%	Yes	pH 2-8
SiO2	5/7/10/20/50/75 $\mu$ m	120 $\text{\AA}$	330m <sup>2</sup> /g	-	-	pH 2-8
SiO2-S	5/10 $\mu$ m	80/100/200 $\text{\AA}$	300m <sup>2</sup> /g	-	-	pH 2-8
SiO2-Bio	5/10 $\mu$ m	300 $\text{\AA}$	120m <sup>2</sup> /g	-	-	pH 2-8

CREDCHROM Preparative HPLC series are versatile media based on the silica-gel substrate for reversed-phase chromatography CredChrom is a long carbon chain C18 bound with the spherical silica-gel substrate which has low metal-ion content ( $<20$  ppm), high specific surface area and good mechanical strength. With CredChrom unique chemical bonding technique, CredChrom has excellent stability of bonding and end-capping. CredChrom meet the high requirements for analysis and preparative applications.


### Advantages:

- Low silanol activity
- Uniform ligand binding
- Low metal content
- Narrow particle size
- Excellent stability


Distribution of Particle Size for CredChrom C18 5um



Repeated Injection tests for CredChrom C18 5um



## National Institute of Standards and Technology (NIST) ERM 870 Test





CREDCHROM C18H can use in alkali environment with high pH CIP (Clean-in-Place) process. The isolation of toluene and ethylbenzene test shows the uniformities of binding ligands on the silica-gel substrate.

# redChrom C18M

## Parameters



Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10um	120Å	330m <sup>2</sup> /g	16%	2-8

## Application




**Nucleotide**  
**Column:** C18M 5µm 4.6×150mm  
 Competitor ODS 5µm 4.6×150mm  
**Mobile Phase:** phosphoric acid buffer / methyl alcohol  
**Flow Rate:** 1ml/min  
**Wavelength:** 254nm  
**Temp.:** 25°C

1 5'-cytidylic acid; 2 5'-uridylic acid;  
 3 5'-guanylic acid; 4 5'-inosinic acid;  
 5 5'-adenylic acid



**Paraben**  
**Column:** EF-C18M 5µm  
 4.6×150mm  
 Competitor ODS 5µm  
 4.6×150mm  
**Mobile Phase:** Water / methyl alcohol  
**Flow Rate:** 1ml/min  
**Wavelength:** 254nm  
**Temp.:** 25°C

1 Methyl ester; 2 Ethyl ester;  
 3 Propyl ester; 4 Butyl ester



**Water-soluble multivitamin**  
**Column:** C18M 5µm  
 4.6×150mm  
**Mobile Phase:** phosphoric acid buffer / acetonitrile  
**Flow Rate:** 1ml/min  
**Wavelength:** 210nm  
**Temp.:** 25°C



1 Pyridoxine;  
 2 VB1;  
 3 Nicotinamide;  
 4 Folic acid;  
 5 VB2

# CredChrom C18H



## Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	120Å	330m <sup>2</sup> /g	20%	2-11

## Application




**Ibuprofen/Benzene ketone**  
**Column:** EF-C18H 5µm 4.6×150mm  
**Competitor** 5µm 4.6×150mm  
**Mobile Phase:** phosphoric acid buffer / acetonitrile  
**Flow Rate:** 2ml/min  
**Wavelength:** 214nm  
**Temp.:** 30°C



**Polar/Nonpolar/Neutral/Alkali Compounds**  
**Column:** EF-C18H 5µm 4.6×250mm  
**Competitor** 5µm 4.6×250mm  
**Mobile Phase:** phosphoric acid buffer / methyl alcohol  
**Flow Rate:** 1ml/min  
**Wavelength:** 254nm  
**Temp.:** 30°C


- 1 Uracil; 2 Butyl p-hydroxybenzoate;
- 3 Propranolol; 4 Di-propyl ortho-phthalate;
- 5 Naphthalene; 6 Acenaphthene;
- 7 Amitriptyline

## The purification of EPA in fish oil



**EPA in fish oil**  
**Column:** C18H 8µm 20×250mm  
**Sample:** 90% EPA material

**Finished sample**  
**Purification:** 99.7%






# CredChrom C8M

## Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	120Å	330m <sup>2</sup> /g	12%	2-8


## Application



**Orlistat**  
**Column:** EP-C8 10µm 10x250mm  
**Mobile Phase:** EtOH solution **Flow Rate:** 4ml/min **Wavelength:** 195nm

**Sample:**  
 Dissolved raw material with methyl alcohol  
**Concentration:** 50-60mg/ml

**Finished sample**  
**Purification:** 99.8%  
**Single impurity** < 0.1%  
**Recovery:** ≥90%





### Insulin

Column: C8 8 $\mu$ m 10 $\times$ 250mm

Time	A	B
0	85%	15%
5min	85%	15%
15min	64%	36%
225min	34%	66%

CredChrom C8	Cycle	Injection	Purification	P1	P1c	P2
	1	100ml	99.76%	0.21%	0.02%	0.01%
		50ml	99.74%	0.22%	0.02%	0.02%
	2	50ml	99.75%	0.22%	0.02%	0.01%
	3	50ml	99.74%	0.22%	0.02%	0.01%
	4	50ml	99.74%	0.22%	0.02%	0.01%
	5	50ml	99.76%	0.21%	0.02%	0.01%
	6	50ml	99.75%	0.22%	0.02%	0.02%
	7	50ml	99.76%	0.21%	0.02%	0.02%
	8	50ml	99.74%	0.22%	0.02%	0.01%
9	50ml	99.74%	0.22%	0.02%	0.02%	

Preparative HPLC Column Size of C18, C8 and Si Phase	Column Size (mm)	Part Number
<b>5um particles</b>		
CredChrom EF-SiO2 5um 120A, pH 2-8	7.8× 250 mm	CCEF-5Si-120-7.8x250
	10× 250 mm	CCEF-5Si-120-10x250
	20× 250 mm	CCEF-5Si-120-20x250
	21.2× 250 mm	CCEF-5Si-120-21.2x250
CredChrom EF-C18M 5um 120A carbon content 16 , pH 2-8	7.8× 250 mm	CCEF-5C18M-120-7.8x250
	10× 250 mm	CCEF-5C18M-120-10x250
	20× 250 mm	CCEF-5C18M-120-20x250
	21.2× 250 mm	CCEF-5C18M-120-21.2x250
CredChrom EF-C18H 5um 120A carbon content 20 , pH 2-11	7.8× 250 mm	CCEF-5C18H-120-7.8x250
	10× 250 mm	CCEF-5C18H-120-10x250
	20× 250 mm	CCEF-5C18H-120-20x250
	21.2× 250 mm	CCEF-5C18H-120-21.2x250
<b>8um particles</b>		
CredChrom EP-C18M 8um 120A carbon content 16 , pH 2-8	7.8× 250 mm	CCEP-8C18M-120-7.8x250
	10× 250 mm	CCEP-8C18M-120-10x250
	20× 250 mm	CCEP-8C18M-120-20x250
	21.2× 250 mm	CCEP-8C18M-120-21.2x250
	30× 250 mm	CCEP-8C18M-120-30x250
	50× 250 mm	CCEP-8C18M-120-50x250
CredChrom EP-C18H 8um 120A carbon content 20%, pH 2-11	7.8× 250 mm	CCEP-8C18H-120-7.8x250
	10× 250 mm	CCEP-8C18H-120-10x250
	20× 250 mm	CCEP-8C18H-120-20x250
	21.2× 250 mm	CCEP-8C18H-120-21.2x250
	30× 250 mm	CCEP-8C18H-120-30x250
	50× 250 mm	CCEP-8C18H-120-50x250
<b>10um particles</b>		
CredChrom EP-SiO2 10um 120A, pH 2-8	7.8× 250 mm	CCEP-10Si-120-7.8x250
	10× 250 mm	CCEP-10Si-120-10x250
	20× 250 mm	CCEP-10Si-120-20x250
	21.2× 250 mm	CCEP-10Si-120-21.2x250
	30× 250 mm	CCEP-10Si-120-30x250
	50× 250 mm	CCEP-10Si-120-50x250
CredChrom EP-C18M 10um 120A carbon content 16 , pH 2-8	7.8× 250 mm	CCEP-10C18M-120-7.8x250
	10× 250 mm	CCEP-10C18M-120-10x250
	20× 250 mm	CCEP-10C18M-120-20x250
	21.2× 250 mm	CCEP-10C18M-120-21.2x250
	30× 250 mm	CCEP-10C18M-120-30x250
	50× 250 mm	CCEP-10C18M-120-50x250
CredChrom EP-C18H 10um 120A carbon content 20 , pH 2-11	7.8× 250 mm	CCEP-10C18H-120-7.8x250
	10× 250 mm	CCEP-10C18H-120-10x250
	20× 250 mm	CCEP-10C18H-120-20x250
	21.2× 250 mm	CCEP-10C18H-120-21.2x250
	30× 250 mm	CCEP-10C18H-120-30x250
	50× 250 mm	CCEP-10C18H-120-50x250